

The Novel Anti-BICD2 Autoantibody Potentially Predicts a Favorable Disease Course in Systemic Sclerosis

Johannes Schulte-Pelkum¹, Petra Budde¹, Hans-Dieter Zucht¹, Daniel Wirtz¹, **Peter Schulz-Knappe¹**, Oliver Distler², Suzana Jordan³, Nicolas Hunzelmann⁴, Matthias Schneider⁵ and Britta Maurer²,
¹Protagen AG, Dortmund, Germany, ²Division of Rheumatology, University Hospital Zurich, Switzerland, ³Department of Rheumatology, University Hospital Zurich, Switzerland, ⁴Department of Dermatology, University of Cologne, Germany, ⁵Rheumatology, Heinrich-Heine-University Düsseldorf, Germany


ACR2016 # 815

Introduction


Systemic sclerosis (SSc) is a systemic autoimmune disease that manifests as progressive fibrosis of the skin and internal organs. SSc is associated with the presence of several autoantibodies to intracellular targets, with the three most important SSc-specific being anti-centromere antibodies, anti-Scl70 antibodies and anti-RNA polymerase III antibodies, which occur in over 50% of SSc patients. Autoantibody specificities are strongly associated with pattern of organ involvement and disease outcome, making autoantibodies an essential tool in the clinical management of SSc. This highlights the need for additional specific and sensitive diagnostic and prognostic biomarkers in SSc. We have recently conducted high-content autoantibody profiling studies of SSc, systemic autoimmune diseases (AID), and healthy controls and found novel SSc-associated autoantibodies. These novel SSc-associated autoantibody biomarker candidates and their diagnostic value were evaluated by testing samples derived from 2 different SSc cohorts.

Clinical annotation data including age, gender, disease duration, MRSS, detailed information on organ involvement making up to a total of either >50 or >100 clinical data points (depending on the site) was available for more than 80% of all SSc serum samples included in this study. All samples were analyzed on the novel Multilisa BICD2 (CE), an ELISA for the semi-quantitative detection of anti-BICD2 antibodies in human serum or plasma.


Within this line of evidence, anti-BICD2 autoantibodies showed a significantly higher titer in patients with moderate skin involvement reflected by low MRSS-stages (p=0.002). Analysis of pulmonary involvement revealed elevated anti-BICD2 titers in the group of patients *not* suffering from pulmonary fibrosis measured by high resolution computer tomography (HCRT, p=0.0001, not shown) or by plain x-ray analysis (p=0.062).




Box-and-whisker plot of anti-BICD2 autoantibody reactivity in SSc, AID, and healthy control serum samples. Elevated reactivity against BICD2 was observed in SSc patients.



ROC analysis of the Multilisa BICD2 tested in two cohorts. SSc sera were identified with a sensitivity of 28.6% and 97.3% specificity (OR=14.56).




Cologne cohort Scatter plot of reactivity of the novel antigen BICD2 compared to staged residual transfer capability lung function in SSc patients. Anti-BICD2 autoantibodies are elevated in the group with low impairment of lung transfer function.



Eustar cohort Scatter plot of reactivity of the novel antigen BICD2 compared to staged residual transfer capability lung function in SSc patients. Anti-BICD2 autoantibodies are elevated in the groups with lowest impairment of lung transfer function.

Target discovery, validation and development: From discovery to IVD assay in three years




The Protagen Biomarker Discovery workflow for the development of novel biomarker targets. SeroTag® autoantibody discovery process is based on a bead-based array screening of recombinant human autoantigens. Assay development is based on ELISA technology. Time from discovery to Product was three years.

Methods


Serum samples were obtained from the biobanks of the Department of Rheumatology, University Hospital Zurich, Switzerland, and Department of Dermatology, University of Cologne, Germany. The analysis included samples collected from patients suffering from SSc (n=301) systemic lupus erythematosus (n=39), Sjögren's syndrome (n=11), rheumatoid arthritis (n=19), idiopathic Inflammatory Myopathy (IIM; n=20) and healthy volunteers (n=99).

Results


We found anti-BICD2 with a prevalence of 28.6% in SSc patients of both cohorts, and only 2.7% prevalence in cohorts with other rheumatic diseases or healthy controls (OR=14.56). Anti-BICD2 autoantibodies were present in a subgroup of SSc patients where skin fibrosis was either from the limited cutaneous subtype (p=0.0002) or restricted to sclerodactyly (p=0.0037).




Scatter plot of anti-BICD2 autoantibody reactivity in the context of modified Rodnan skin score. High reactivity against BICD2 was observed in SSc patients suffering only from mild impairment of the skin as reflected by stages 0 and 1, statistical significance calculated by Kruskal Wallis test.



Scatter plot of anti-BICD2 autoantibody reactivity in the context of skin involvement. High reactivity against BICD2 was observed in SSc patients suffering from limited skin involvement or from sclerodactyly.



Elevated anti-BICD2 reactivity is found in patients with negative diagnosis of pulmonary fibrosis (by X-Ray)



Box plot of disease duration. Anti-BICD2 reactivity was found higher in the group with disease duration of more than 3 years.

Conclusion

In this study, we were able to further confirm the diagnostic value and high specificity of the newly discovered BICD2 autoantigen. Anti-BICD2 autoantibodies were found to be elevated in patients suffering from limited SSc, and anti-BICD2 reactivity was found to be related to clinical observations of a moderate course of disease. The observations of anti-BICD2 being found in patients with limited cutaneous involvement, low MRSS, absence of pulmonary fibrosis and an only moderate impairment of lung function point towards a moderate course of disease in anti-BICD2 positive SSc patients.