

# Development of a qualitative ELISA for the detection of anti-BICD2 autoantibodies in systemic sclerosis

Daniel Wirtz<sup>1</sup>, Johannes Schulte-Pelkum<sup>1</sup>, Petra Budde<sup>1</sup>, Peter Schulz-Knappe<sup>1</sup>, Hans-Dieter Zucht<sup>1</sup>, Suzana Jordan<sup>2</sup>, Oliver Distler<sup>2</sup>, Britta Maurer<sup>2</sup>, Nicolas Hunzelmann<sup>3</sup>, Karsten Conrad<sup>4</sup>, Matthias Schneider<sup>5</sup>

<sup>1</sup>Protagen AG, Dortmund, Germany; <sup>2</sup>Division of Rheumatology, University Hospital Zurich, Switzerland; <sup>3</sup>Universitätsklinikum Köln, Klinik für Dermatologie und Venerologie, Köln, Germany; <sup>4</sup>Carl Gustav Carus Universität, Dresden, Germany; <sup>5</sup>Heinrich Heine Universität, Poliklinik für Rheumatologie, Düsseldorf, Germany

## Introduction

Systemic sclerosis (SSc) is a systemic autoimmune disease that manifests as progressive fibrosis of the skin and internal organs<sup>1,2</sup>. We recently conducted autoantibody profiling studies in SSc using the Protagen SeroTag<sup>®</sup> process and found novel SSc-associated autoantibodies with the potential to improve the diagnosis of SSc and other systemic autoimmune diseases (SARD). Beyond the classical markers such as anti-centromere and anti-Scl70 antibodies, which are generally exclusive of the other, novel autoantibodies were identified. Here, we report the development of the development of a qualitative ELISA for the detection of human BICD2 autoantibodies. BICD2 (Protein bicaudal D homolog 2) is an evolutionarily conserved motor adaptor protein that is involved in dynein-mediated transport by linking the dynein motor complex to various cargoes.

## Methods

We developed an ELISA using a highly purified recombinant fragment of human BICD2, the Protagen<sup>®</sup> MULTILISA<sup>®</sup> BICD2-ELISA.

The assay principle is described in Fig. 1.

Two lots of kit-material were independently manufactured using two different lots of purified BICD2 protein. Quality control was performed with a quality control panel of 8 samples.

The MULTILISA<sup>®</sup> BICD2-ELISA was calibrated using receiver operating characteristics (ROC) analysis to 95% specificity.




Figure 1: Test procedure of the Protagen<sup>®</sup> MULTILISA<sup>®</sup> BICD2-ELISA

## Diagnostic performance

A mixed cohort consisting of SSc, SARD and healthy donor samples (n=394) were tested with the Protagen<sup>®</sup> MULTILISA<sup>®</sup> BICD2-ELISA. Data were analyzed using ROC analysis (Fig. 2).




Figure 2: ROC analysis of 394 samples of a mixed cohort tested with Protagen<sup>®</sup> MULTILISA<sup>®</sup> BICD2-ELISA. The ROC analysis show 31% sensitivity at 95% specificity and 22% sensitivity at 99% specificity (Table 2). The cutoff was determined at 0.650 OD and fixed by a ready-to-use calibrator.

## Linearity analysis

Four positive samples were stepwise diluted with a negative sample and measured using the Protagen<sup>®</sup> MULTILISA<sup>®</sup> BICD2-ELISA. The calculated Index (ratio from OD<sub>sample</sub> and OD<sub>calibrator</sub>) was plotted over sample amount (Fig. 3).

The samples show linear dilution behavior in a range between I=0.2 and I=4. This covers an OD Range between 0.150 OD and 2.600 OD.




Figure 3: Dilution series of 4 positive sample with a negative sample. Plotted are the calculated Index of the sample against the amount of sample in double logarithmic scale

## Added value

173 SSc-samples were analyzed with one lot of the Protagen<sup>®</sup> MULTILISA<sup>®</sup> BICD2 ELISA. SSc patients tested positive for anti-BICD2, anti-Scl-70 and anti-CENPB autoantibodies are shown in a Venn-Diagram (Fig. 3).

The prevalence of a-Scl70, a-CENPB and a-BICD2 autoantibodies were comparable in this cohort. The anti-BICD2 assay identifies 7 additional patients, who were tested negative for anti-Scl70 and anti-CENPB.




Figure 3: VENN Diagram showing the number of SSc sera with positive results for BICD2, Scl70 and CENPB

## Verification data

Two lots of material were used to verify the assay and to determine the lot-to-lot consistency, intra- and inter-assay-variability, robustness against interference substances and the accelerated stability.

### Lot-to-lot consistency

Repeatability n=72	Within Lot			Lot-to-Lot		Cal-to-Cal	
	Slope	95% CI	Intercept	95% CI	95% CI	95% CI	95% CI
	1.09	1.04 to 1.11	0.006	0.95 to 1.03	0.007	-1.59*10 <sup>-17</sup> to 2.1*10 <sup>-17</sup>	0.958 to 0.958

### Intra-/inter-assay variance

Intra-assay n=30						Inter-assay n=150									
Unit = Index (Sample/Calibrator)						Unit = Index (Sample/calibrator)									
Sample	Mean	SD	95% CI	CV		Mean	SD	95% CI	CV		Mean	SD	95% CI	CV	
Sample 1 (positive)	1.6	0.1	0.0	0.2	5%	Sample 1 (positive)	1.6	0.1	0.1	0.1	6%				
Sample 2 (negative)	0.4	0.0	0.0	0.0	2%	Sample 2 (negative)	0.4	0.0	0.0	0.0	7%				
Sample 3 (negative)	0.5	0.0	0.0	0.0	4%	Sample 3 (negative)	0.5	0.1	0.0	0.1	11%				
Sample 4 (Intermediate)	1.1	0.1	0.1	0.2	8%	Sample 4 (Intermediate)	1.2	0.2	0.2	0.3	16%				
Sample 5 (high positive)	3.1	0.1	0.0	0.2	2%	Sample 5 (high positive)	3.2	0.1	0.1	0.2	4%				

### Robustness against interference substances

Interference substance (n=108)	Result
Hemoglobin	Up to 300mg/dl
Bilirubin	Up to 18 mg/dl
Intralipid	Up to 1000mg/ml
Rheumatoid factor	Up to 500 IU/ml

### Accelerated stability at 37°C

Residual activity [%]	Plate at 37°C [days]	Calibrator at 37°C [days]	Conjugate at 37°C [days]
90%	19	52	20
85%	27	68	30
80%	35	83	41

## Assay portability studies

Test portability was verified using bead based assay, and Western blot / Line immunoassay studies. The preliminary results show high correlation of all results to ELISA (data not shown).

## Conclusion

We developed a qualitative ELISA for detecting anti-BICD2 autoantibodies, which occur in 30% of SSc patients. The assay was validated by analyzing SSc samples and SARD control samples. The availability of a validated, qualitative BICD2-ELISA may improve the specific detection of SSc and reduce the diagnostic gap.

Verified RuO-Kits of the Protagen<sup>®</sup> MULTILISA<sup>®</sup> BICD2-ELISA will be available in November 2015.

### References

- Tan EM: Autoantibodies to nuclear antigens (ANA): their immunobiology and medicine. *Adv Immunol* 1982, 33:167-240.
- Hudson M, Fritzler MJ: Diagnostic criteria of systemic sclerosis. *J Autoimmun* 2014, 48-49:38-41.
- Meier FMP, Frommer KW, Dinser R, Walker U a., Czirjak L, Denton CP, Allanore Y, Distler O, Riemekasten G, Valentini G, Muller-Ladner U: Update on the profile of the EUSTAR cohort: an analysis of the EULAR Scleroderma Trials and Research group database. *Ann Rheum Dis* 2012, 71:1355-1360.
- Müller-Ladner U, Tyndall A, Czirjak L, Denton C, Matucci-Cerinic M: Ten years EULAR Scleroderma Research and Trials (EUSTAR): what has been achieved? *Ann Rheum Dis* 2014, 73:324-7.